



Imaging and Flow Cytometry Core

NovoCyte Advanteon BVR Standard Operation Protocol Basic Operation

1. NovoExpress Software Log In

Log into NovoExpress software with your own login name and password. Make sure *Auto Login* is unchecked.

**Please contact Faculty Core Facility Staff to establish a new user account.*

Group: Organization
Username: administrator
Password: [empty]
 Auto Login
Login Quit

*Username *Ocleaning* do not have password. Leave password entry box empty and click *Login*.

Ocleaning account is for cleaning purpose only, **DO NOT USE THIS ACCOUNT TO PERFORM EXPERIMENT.**




Imaging and Flow Cytometry Core

2. Compensation (Perform when needed)

Step 1. Select appropriate **plate type** in the *Plate Manager*.

Step 2. Set-up Compensation Controls

a. In the *Home* tab of the **Menu Bar**, click the *Auto Compensation* button.



b. Select Compensation on: *Height*, Parameter for calculation: *Median* and **check the boxes of channels involved**. Then click **OK**



Imaging and Flow Cytometry Core

Import Samples from FCS files

Compensation on:

Area
 Height

Parameter for calculation:

Mean
 Median

Compensation Channels:

Unstained All [Reset All Photodetector Gain](#)

<input checked="" type="checkbox"/> FSC	FSC	<input checked="" type="checkbox"/> R660	APC	<input checked="" type="checkbox"/> V695	Qdot 705
<input checked="" type="checkbox"/> SSC	SSC	<input checked="" type="checkbox"/> R695	Alexa Fluor 680	<input checked="" type="checkbox"/> V725	BV711
<input checked="" type="checkbox"/> B530	FITC	<input checked="" type="checkbox"/> R725	Alexa Fluor 700	<input checked="" type="checkbox"/> V780	Qdot 800
<input checked="" type="checkbox"/> B586	EYFP	<input checked="" type="checkbox"/> R780	APC-Cy7	<input checked="" type="checkbox"/> Y586	PE
<input checked="" type="checkbox"/> B615	PI	<input checked="" type="checkbox"/> V445	Pacific Blue	<input checked="" type="checkbox"/> Y615	PE-Texas Red
<input checked="" type="checkbox"/> B660	PerCP	<input checked="" type="checkbox"/> V530	AmCyan	<input checked="" type="checkbox"/> Y660	PE-Cy5
<input checked="" type="checkbox"/> B695	PerCP-Cy5.5	<input checked="" type="checkbox"/> V586	Pacific Orange	<input checked="" type="checkbox"/> Y695	PE-Cy5.5
<input checked="" type="checkbox"/> B725	PerCP-eFluor 7	<input checked="" type="checkbox"/> V615	Qdot 605	<input checked="" type="checkbox"/> Y725	PE-Alexa Fluor
<input checked="" type="checkbox"/> B780	PE-Cy7 (B)	<input checked="" type="checkbox"/> V660	Qdot 655	<input checked="" type="checkbox"/> Y780	PE-Cy7 (Y)


OK Cancel

- c. Compensation Control Specimen is created in the Experiment Manager panel with corresponding empty control samples of specific position of the tube rack or plates (e.g. A3: B530 FITC, A3 is the position of the rack or plate and refers to FITC single stain controls).

The compensation controls tubes should be placed in the rack or plate according to the positions given (i.e. Put FITC single stain tubes in A3 position).



Imaging and Flow Cytometry Core




You may change the position by right-click of the sample name and *Rename* the control (e.g. A3: B530 FITC can be renamed to B3: B530 FITC, FITC single stain tube position now changes from A3 to B3).

Imaging and Flow Cytometry Core

Step 3. Compensation Control Acquisition

- a. Click **Run Plate** on the Cytometer Control Panel.




- b. Select wells of compensation controls. Click **Run** and then **OK** to proceed.
- c. After all controls have been acquired, the compensation matrix is calculated automatically.

Step 4. Apply Compensation Matrix to Experiment Sample

- a. **Drag** the *Compensation* node under the *Compensation Specimen* and **Drop over** the desired sample.



Imaging and Flow Cytometry Core





CPOS - Imaging and Flow Cytometry Core



Imaging and Flow Cytometry Core



- b. To fine tune the Compensation, click on the plot you want to adjust and click the *Quick Compensation* button In the *Home* tab of the Menu Bar OR the *quick compensation icon* in the tool bar.




Scrollbars appear on any two parameters plots with fluorescent parameters opened on the workspace. Quickly adjust compensation by **dragging the scrollbar**.



Imaging and Flow Cytometry Core



c. To view or adjust the compensation matrix, lick the **Compensation Matrix** button In the **Home** tab of the **Menu Bar** and the Compensation Matrix window will show.





Imaging and Flow Cytometry Core

The screenshot shows a software window titled "Compensation for T-lym". It contains a table with two tabs: "Compensation Matrix" (selected) and "Spillover Matrix". The table lists various fluorescent markers and their compensation values. The "OK" button at the bottom right is highlighted with a red box.

Source\Target	CD3 FITC	CD8 PE	CD45 PerCP	PE-Cy7	CD4 APC	APC-Cy7
CD3 FITC	-100.3801	8.9548	-2.1786	0	0.4603	0
CD8 PE	4.2611	-100.3911	27.5831	0	-5.84	0
CD45 PerCP	0.0185	0.0386	-100.6337	0	21.6362	0
PE-Cy7	0	0	0	-100	0	0
CD4 APC	-0.0005	-0.0011	2.8983	0	-100.6231	0
APC-Cy7	0	0	0	0	0	-100

Below the table, there is a slider for "- % CD3 FITC" and several buttons: "Preview", "Clear", "Restore", "OK", and "Cancel". The "OK" button is highlighted with a red box.

To adjust, **check *Preview*** box and adjust the corresponding value. The corresponding plots will refresh with updated value real time. Adjust until satisfied. Then click ***OK*** to apply.

Click *Restore* to restore the Auto-compensation matrix value.

CPOS - Imaging and Flow Cytometry Core




Imaging and Flow Cytometry Core

3. Sample acquisition with NovoSampler Q

Step 1. Create experiment samples from the Plate Manager

a. Select appropriate Plate type. Choose **40-tube rack** for 5-mL flow tubes




b. Highlight the position with samples on the plate by holding left Click and Drag AND/OR hold Ctrl and left-click to multi-select specific wells. Black square indicates selected well.



Imaging and Flow Cytometry Core

	1	2	3	4	5	6	7	8	9	10	11	12
A		1	2	2	2	2						
B	3	3	4	4	5	5	6	6	7	7	8	8
C	9	9	10	10	11	11	12	12	13	13	14	14
D	15	15	16	16	17	17	18	18	19	19	20	20
E	21	21	22	22	23	23	24	24	25	25	26	26
F												

c. Click the *New Sample(s)* button to create a new sample of *Specimen 1*.



d. Repeat step *1b* and *1c* to create new sample of *Specimen 2* if needed.


e. Check *Absolute count* if absolute counting is required.



Imaging and Flow Cytometry Core

*Dead volume will increase from 10 μ L to 30 μ L with Absolute count checked.

- f. Double click *Sample 1* on the Experiment Manager until the red arrow is pointing to Sample 1.



Step 2. Select Channels

- a. Click on the “A” and “H” of the parameters panel in Cytometer setting to Select OR Unselect ALL.



Imaging and Flow Cytometry Core

Cytometer Setting

Parameters:

Pa...	Alias	Gain	A	H
FSC	FSC	364	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SSC	SSC	364	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B530	FITC	462	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B586	EYFP	542	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B615	PI	520	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B660	PerCP	525	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B695	PerCP-Cy5.5	678	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B725	PerCP-eFluor 710	283	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B780	PE-Cy7 (B)	354	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Stop Condition:

12,000 Events on

0 Min Sec

50 μ L

- b. *Check the box of A or H* of the interested channels to select. Please always check A for FSC (H is checked by default).

Cytometer Setting

Parameters:


Pa...	Alias	Gain	A	H
FSC	FSC	364	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SSC	SSC	364	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B530	FITC	462	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B586	EYFP	542	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B615	PI	520	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B660	PerCP	525	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B695	PerCP-Cy5.5	678	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B725	PerCP-eFluor 710	283	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B780	PE-Cy7 (B)	354	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

Step 3. Conditions Setup

Imaging and Flow Cytometry Core

- a. Set up the data recording stop conditions by **checking the box next to the condition** *Events* and/or *Time* and/or *Volume*. Acquisition will stop when ANY one of the selected condition(s) is fulfilled.

*Volume is compulsorily selected.



Events

Time

Volume


Range of each conditions:

Events	1 – 10,000,000
Time	0-60 min; 0-59 Sec
Volume	5 – 5000 µL

* Volume may be limited to the Plate format. Please refer to the Appendix.

- b. Select flow rate by click the radio button of *Slow* (14 µL/min), *Medium* (35 µL/min), and *Fast* (66 µL/min) OR use the slider to adjust the flow rate from 5 ~120 µL/min.

* Current sample's flow rate and the corresponding core diameter are shown in the bottom of the panel.





Imaging and Flow Cytometry Core

- c. Set the appropriate threshold by select the appropriate parameters and type in the appropriate number on the *Threshold* panel.

Threshold: Adjust on Plot

FSC-H larger than 100,000

- larger than 10

Storage Gate Ungated

Suggested Threshold on Different cell type:

Cell Types	FSC-H Threshold
Cell lines, larger than 20 μm in cell diameter	300,000~1,000,000
Cells lines, smaller than 20 μm in cell diameter	100,000~300,000
Fixed or un-fixed, freshly isolated cells (leukocytes, spleen cells, thymocytes)	50,000~200,000
Platelets	5,000~10,000
Bacteria	1,000~10,000

- d. Setup mixing and rinsing conditions under Plate Manager.

*For 96-well plate, the default is 1000 rpm.



Imaging and Flow Cytometry Core

Mix Every Well(s)

Rinse Every Well(s)

Mixing Parameters

Speed rpm Acceleration s


Duration s

CPOS - Imaging and Flow Cytometry Core

Imaging and Flow Cytometry Core

Step 4. Draw Plots

a. Click the icon of the interested **plot type** above the workspace



Plot type	Number of parameters	Description
Dot plot	2	The intensities of two parameters are represented by the coordinates of an event (one dot) on the plot.
Density plot	2	The intensities of two parameters are represented by the coordinates of an event (one dot) on the plot with colour-coded density display.
Contour plot	2	The intensities of two parameters are represented by the coordinates on the plot with contour line to show density.



Imaging and Flow Cytometry Core


Histogram plot	1 (x axis only)	The intensity of a parameter is represented along the x-axis, and the number of events at each intensity value is represented along the y-axis.
----------------	--------------------	---

CPOS - Imaging and Flow Cytometry Core



Imaging and Flow Cytometry Core


- b. Create the following plots with the following sequence.
- FSC-H VS SSC-H (Mother population of interest) >
 - FSC-H VS FSC-A (Single Cell Gate) >
 - Live-Dead VS SSC-A (if applicable) >
 - Fluorescence Plots (if applicable)





Imaging and Flow Cytometry Core


- c. To change the parameters of a plot, **mouse over the axis label** and **right-click** to open the drop-down menu of parameters list. Select the parameter of interest.



- d. **Right-click** within a plot to change the plot type if needed.




Imaging and Flow Cytometry Core




Imaging and Flow Cytometry Core

- e. To copy all the settings and plots to other samples in Specimen 1, **drag Sample 1** and drop over *Specimen 1* on the Experiment Manager.



- f. Click *Paste to All*.






Imaging and Flow Cytometry Core

Step 5. Save Experiment

- a. Click **File** on the Menu bar.



- b. Click *Save As*.



- c. Save the experiment (.ncf) in the folder below.

Computer> Experiment Data (D:)> Users> Department> YOUR FOLDER

*Default location: Computer> D: > Novoexpress Data. **Please do NOT save as default**

Click *Save*.

Imaging and Flow Cytometry Core

- Step 6. Load tube rack/ plate
- Lift the cover of the NovoSampler Q.



- Place the tube rack with your sample tubes or plate on the orbital Shaker with A1 position on the top left-hand corner. Make sure the rack or plate is placed within the 4 metal poles.




- Put down the cover.



Imaging and Flow Cytometry Core

Step 7. Sample Acquisition


- a. Click **Mix Sample** to perform orbital shaking.




- b. To run **SINGLE** well/ tube, **double-click** on the interested well on Plate Manager and the selected position highlight in red.
*A2 is selected in the picture.



Imaging and Flow Cytometry Core




- c. Click **Run Single Well** on the Cytometer Control Panel.





Imaging and Flow Cytometry Core

- d. To run multiple tubes / wells automatically, click **Run Plate** on the Cytometer Control Panel.




- e. Double check if the plate type is correct as it state. Click **Run** to continue, or else click **Cancel** and correct the plate type in Plate Manager panel.



Imaging and Flow Cytometry Core

- f. **Select the tubes or wells** you would like to be acquired on the Plate View. Selected wells highlight in Blue. Then click **Run**.




- g. Click **OK** to continue.

Imaging and Flow Cytometry Core

4. Data Analysis during acquisition


Step 1. Set the appropriate display range of the plot.

- Select the FSC-H Vs SSC-H plot (The colour of the header of the plot will be darker). Click **Auto range button** to optimize the data display range.




- To fine tune the data display range, click **zoom in / zoom out buttons**.

Imaging and Flow Cytometry Core




- c. **Drag** on the interested region on the plot if you click *zoom in*.



- d. **Click within a plot** if you click *zoom out*. The range increases by 20% of the current range. Click repeatedly until the desired range is reached.
- e. To change the scale of parameters, **right click on the coordinate label** to open and select the axis scaling (i.e. Linear, Log or Biexponential). Click **Setting** for more options.

Imaging and Flow Cytometry Core



- f. If you cannot achieve a desirable range by using the plot range tools, **adjust the Gain** of the corresponding channels in Cytometer Control - Parameters.

To adjust photodetector gain of one parameter, **double click the current Gain number** of the specified parameter, the photodetector gain adjustment slider will show. **Drag the slider bar or directly enter the value** to change the photodetector gain.

***Gain can only be adjusted during acquisition.**


Parameters:			
Pa...	Alias	Gain	A H
FSC	FSC	364	+ /
SSC	SSC	364	/
B530	FITC	462	/
B586	EYFP	643	/
B615	PI	520	/
B660	PerCP	525	- / Reset /










Imaging and Flow Cytometry Core

Step 2. Gating

- Draw Gates to gate out the target population on the FSC-H VS SSC-H plot with gating tools.



rectangular gates , elliptical gates , polygonal gates , quadrant gates , logic gates ,
range gates , and bi-range gates 


Dot Plot, Density Plot, Contour Plot – All gates suitable

Histogram Plot – Range / Bi-range gate suitable


Imaging and Flow Cytometry Core

- b. To create *rectangular/ elliptical/ range/ bi-range gate*, click the corresponding icon and drag in the plot to enclose the target population within the shape. Release the mouse button to create the gate.

To create *polygonal gate*, click the corresponding icon and left click in the plot to create the first vertex of the polygon. Click in a new location to create the second vertex of the polygon. Continue moving around the target population and creating vertices until the target population is enclosed. On the last vertex, double-click to complete the polygon and create the gate.




To create *quadrant gate*, click the corresponding icon and Click in the plot to create the center of the quadrants and create the gate. As shown below, the center, endpoints, and lines of the quadrant gate can be moved to enclose the correct populations.



Imaging and Flow Cytometry Core

- c. To create gate subpopulation, **right-click at the plot header** of a plot to display a drop-down menu. Select the mother gate and create a new gate for your target.





- d. The Gate Manager panel displays all gates of the active sample in list mode or tree mode. It provides user interface to modify gate name, color and color precedence and also shows gate hierarchy and gate statistics.

Gate	Color	Count	% Parent	X	Y	Mean X	Mean Y	CVX	CVY
All	8	10,035							
GR	3	3,160	31.49%	CD45 PerCP-H	SSC-H	3,207	1,017,885	29.33%	15.15%
MO	2	560	5.58%	CD45 PerCP-H	SSC-H	5,727	403,460	22.95%	15.38%
LY	1	4,945	49.28%	CD45 PerCP-H	SSC-H	14,276	109,626	19.14%	27.48%
CD3-CD4+	4	7	0.14%	CD3 FITC-H	CD4 APC-H	975	27,054	56.15%	48.79%
CD3+CD4+	5	1,660	33.58%	CD3 FITC-H	CD4 APC-H	7,125	40,802	34.06%	17.69%
CD3-CD4-	6	2,044	41.34%	CD3 FITC-H	CD4 APC-H	481	5,402	45.91%	22.28%
CD3+CD4-	7	1,233	24.94%	CD3 FITC-H	CD4 APC-H	6,113	5,240	38.80%	24.66%

Imaging and Flow Cytometry Core

Step 3. Statistics


- a. To edit statistics, click the button on the lower right corner of a plot to expand the plot and display the statistics chart first.



- b. Right-click within the chart** and select the parameters to hide or display.




Imaging and Flow Cytometry Core



5. Data Export

Step 1. Export FCS file

- a. Select the sample, specimen, group, or experiment file node with data to be exported in the Experiment Manager Panel. **Right-click** the node and select **Export** → **Export to FCS Files ...** The Export Events window will open.



- b. Choose “All” for the *Gate* option. Click “...” button next to entry box of *Path*. Select your saving destination in Experiment Data Drive (F:/) **Experiment Data (D:)> user > Department> Your NAME**



Imaging and Flow Cytometry Core

Select “*FSC3.0*” for Format and Click *OK*.

- Data files (including the experiment files) can be stored in D drive for **3 months ONLY**. While experiment template files may be stored for a longer duration.
- Data on the computer (D drive) will be removed regularly without prior notice.


CPOS - Imaging and Flow Cytometry Core




Imaging and Flow Cytometry Core

Step 2. Export PDF file (optional)

- a. To Export Plots and Statistics to a PDF, **double-click the Report node** in the Experiment Manager panel and Report Window will popup.




- b. Click the *PDF* button of the tool bar.



Select your saving destination in Experiment Data Drive (D:/)
D:/user/Your department/Your NAME

Imaging and Flow Cytometry Core


- c. Click *Save*  buttons on the top left-hand corner of the window when you finish your experiment

Step 3. Export experiment template file (optional)

- a. To export experiment template, select the experiment file node to be exported in the Experiment Manager Panel. **Right-click** the node and select **Export** → **Export as Experiment Template**.

6. System Cleaning (You may use the “Ocleaning” account to perform)

- a. Place tubes of at least 1ml of cleaning solution 1, 2 and 3 and put them in A1-A3 of the 40-tube rack respectively.




- b. Click *File* on the Menu bar.




- c. Click *New* > *New from Experiment File*.

Imaging and Flow Cytometry Core



- d. Select *Desktop > cleaning.ncf* . Then click **OK**
- e. Click *Run Plate* on the Cytometer Control Panel.



- f. Select all wells. Then click *Run*. (Click **OK** to continue).
- g. Select *Desktop > cleaning.ncf*. Click **Save** and **Yes** to overwrite.



**HKU
Med** LKS Faculty of Medicine
Centre for PanorOmic Sciences
香港大學泛組學科研中心

Imaging and Flow Cytometry Core


CPOS - Imaging and Flow Cytometry Core




Imaging and Flow Cytometry Core

7. Re-use Experiment as template

a. Click *File* on the Menu bar.



b. Click *New* > *New from Experiment File*.



c. Select your target experiment file (.ncf) . Then click *OK*


d. Click *File*> *Save As* to save the new experiment.




Imaging and Flow Cytometry Core

8. NovoExpress Software Log out

a. Click *File* on the Menu bar.



b. Click *Logout*



Imaging and Flow Cytometry Core

APPENDIX

Minimum Sample Volume Requirements for the NovoSampler Q is shown in table 3-1 below

Parameter	Plate Type	Standard Mode/HT Mode	Custom Mode	
			With Absolute Count Checked	With Absolute Count Unchecked
Minimum Sample Volume (μL)	12 \times 75 mm tube (with ACEA 40 tube rack)	20.5	40.5	20.5
	24-well plates	285	305	285
	48-well plates	110	130	110
	96-well plates (flat-bottom)	41.5	61.5	41.5
	96-well plates (V-bottom)	15.6	35.6	15.6
	96-well plates (U-bottom)	15	35	15
	384-well plates	18	38	18

Table 3-1 Minimum Sample Volume Requirement for Each Plate and Tube Type



The recommended sample concentration range is 1×10^6 to 5×10^6 cells/mL.